

Description

The set consists of:

Primary stain, crystal violet, modified

Hucker: Crystal violet 10%, Ethanol, 95%, Ammonium oxalate ((NH₄)₂C₂O₄•H₂O) 1%

Presentation: As a ready-to-use solution or as solution A + solution B to be mixed.

Stain intensifier: Lugol's iodine 0.33%, Iodine (I₂) 0.33%, Potassium iodide (KI) 0.67% ,

Solvent: Distilled water

Decolourizer: Ethanol 95%

Counterstain: Safranine, 0.25%,
Solvent: Ethanol 10%


Grade: For IVD (in-vitro diagnosis)

Supplied as a set of 4 opaquer bottles of each 250ml

Educational Lab Equipments,
#449, HSIIDC, Industrial Area, Saha,
Haryana

Direct Contact Details 

+91-98173-19615  sales@educational-equipments.com

 www.educational-equipments.com