

## Description

---

**The set consists of:**

**Primary stain, crystal violet, modified**

**Hucker:** Crystal violet 10%, Ethanol, 95%, Ammonium oxalate ((NH<sub>4</sub>)<sub>2</sub>C<sub>2</sub>O<sub>4</sub>•H<sub>2</sub>O) 1%

**Presentation:** As a ready-to-use solution or as solution A + solution B to be mixed.

**Stain intensifier:** Lugol's iodine 0.33%, Iodine (I<sub>2</sub>) 0.33%, Potassium iodide (KI) 0.67% ,

**Solvent:** Distilled water

**Decolourizer:** Ethanol 95%


**Counterstain:** Safranine, 0.25%,  
Solvent: Ethanol 10%

Grade: For IVD (in-vitro diagnosis)


Supplied as a set of 4 opaquer bottles of each 250ml

---

Educational Lab Equipments,  
#449, HSIIDC, Industrial Area, Saha,  
Haryana

Direct Contact Details 

+91-98173-19615  sales@educational-equipments.com

 [www.educational-equipments.com](http://www.educational-equipments.com)